



# Changing the wax traps




Changing the Wax Traps	Latest review date: Apr-2024	Version number: 1.4
Status: Approved	Next review: Oct-2024	Page 1 of 4


# Changing the wax traps




1<sup>st</sup> remove the dome



Underneath you will find a white circular grommet this is the wax trap




Take a new wax trap from the pack. On the empty end of the applicator stick is a barbed protrusion. Press this into the wax trap you want to remove.



Remove the old wax trap.


Changing the Wax Traps	Latest review date: Apr-2024	Version number: 1.4
Status: Approved	Next review: Oct-2024	Page 2 of 4



Now flip the applicator around and press the new wax trap in to replace the old one.



Remove the applicator and the new wax trap will stay in place.



Finally replace the dome.

The wax traps should only be replaced if they are blocked.

**Always make sure the dome is safely secured before you put the hearing aid into your ear.**

Changing the Wax Traps	Latest review date: Apr-2024	Version number: 1.4
Status: Approved	Next review: Oct-2024	Page 3 of 4

# Useful information

## Repair Clinics

We offer walk-in (open) repair clinics at two of our Exeter sites, and appointment only (bookable) repair clinics at several of our localities. This is for servicing of hearing aids and simple repairs (e.g. ear mould re-tubing, ear mould modifications and replacement aids).

	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday
<b>Wonford (RD&amp;E)</b>	Open	Open	Open	Open	Open	Bookable
<b>CHC Queen St. Exeter</b>	Open	Open	Open	Open	Open	Open
<b>Budleigh</b>					Bookable	
<b>Crediton 2nd &amp; 4th Friday</b>					Bookable	
<b>Cullompton 2nd Friday</b>					Bookable	
<b>Exmouth</b>			Bookable			
<b>Honiton</b>					Bookable	
<b>Okehampton</b>				Bookable		
<b>Seaton</b>		Bookable				
<b>Sidmouth</b>				Bookable		
<b>Tiverton</b>	Bookable					

**If your personal details change, please inform us directly.**

**Address:** Chime Social Enterprise, Audiology Department,  
Royal Devon & Exeter Hospital (Wonford),  
Barrack Road, Exeter, Devon, EX2 5DW

**Tel:** 01392 402223

**E-Mail:** [cse.audiology@nhs.net](mailto:cse.audiology@nhs.net)

[www.chimehealth.co.uk](http://www.chimehealth.co.uk)

Chime provides NHS Audiology Services for NHS Devon - Mid, East and Exeter areas. From 1st May 2011 the existing audiology department including staff and equipment transferred to Chime - the new Social Enterprise Company. Chime is located at the RD&E Foundation Trust (but is separately managed).

Changing the Wax Traps	Latest review date: Apr-2024	Version number: 1.4
Status: Approved	Next review: Oct-2024	Page 4 of 4